



A β detection, immunoblot

HEK293 APP^{swe} - Huang et al., J Neurosci 2016

Detection of A β in conditioned media by immunoblot requires concentration (using a microcon system, or by TCA) – which requires conditioning in FBS-free media (FBS will interfere with protein concentration methods). As A β monomers are small (<10kDa), precipitates or concentrated media are resolved on 10-20% Tricine gradient gels. Other secreted or shed extracellular factors (sAPP α , sSORLA) are more abundant, and may not require concentration steps for detection by immunoblot.

Steps:

- 1) Wash cells twice with plain DMEM (no FBS) for 4-16h before harvest (incubate in 1ml DMEM ea. 6-well plate; typically media is conditioned overnight).
- 2) TCA precipitate proteins in the conditioned media; lyse cells and immunoblot lysates for proteins of interest.
- 3) Dissolve pellet with Tricine loading buffer.
- 4) Run 10-20% Tricine gel, use Tricine running buffer.
- 5) Transfer (immunoblot) onto nitrocellulose or PVDF.
- 6) Incubate the membrane in 0.2% Glutaraldehyde for 30mins.
- 7) Block membranes with 5% milk for 1h.
- 8) Probe for A β using the 6E10 antibody.

TCA precipitation

Steps

- 1) Add 2% Na-deoxycholate to 0.02% (100x) to the FBS-free conditioned media
- 2) Mix at RT (room temperature) for 15 mins.
- 3) Add 100% TCA to a final concentration of 10%, incubate rocking 1h at RT.
- 4) Spin 4°C for 10 mins. (max speed), remove supernatant and dry.
- 5) Add 200 μ ls of cold acetone to TCA pellet, incubate on ice 15 mins.
- 6) Spin 4°C for 10 mins at max speed.
- 7) Remove supernatant and dry tubes, resuspend in Tricine running buffer.

Note: for larger proteins, standard SDS-PAGE Tris/glycine gels can be used