



Transformation of *E. coli* by CaCl_2

This protocol is used to generate competent cells (XL1, DH5 α) (**part A**) which are then frozen in aliquots at -80°C . 50 μl s is generally sufficient for a single transformation (**part B**) from a ligation reaction.

A. Competent cells:

1. Grow 1 liter of cells to $\text{OD}_{600} = 0.4$ (2x 10ml seeds – 37deg).
2. Cool on ice for 10 mins, spin 20 mins. at 6000 rpm.
3. Resuspend in 500mls of ice cold 50mM CaCl_2 /10mM Tris-HCl pH 7.4 (sterile filtered).
4. Incubate on ice for 1h.
5. Spin, resuspend in 5mls of 50mM CaCl_2 , 10mM Tris-HCl pH 7.4, 20% glycerol, sterile filtered.
6. Aliquot (150 μl s) and freeze at -80°C .

B. Transformation:

1. Thaw cells on ice.
2. Add DNA and 100 μl s TCM, 100 μl s of cells on ice, incubate 30mins.
3. Heat shock at 42°C for 0.5 min, ice for 1.5 mins.
4. Add 2mls of LB, incubate 37°C for 1h.
5. Spin down 3k for 15 mins., plate on LB plus antibiotic.

TCM: 50mM CaCl_2 , 30mM MgCl_2 , 10mM Tris-HCl pH 7.4